

## B10-13, block

### DESCRIPTION

Ø13 mm x 10 sockets, flat bottom, depth 30 mm

### SPECIFICATIONS

|          |       |
|----------|-------|
| Depth    | 30 mm |
| Capacity | 10    |
| Diameter | 13 mm |



### CAT. NUMBER

BS-010418-LK

B10-13



CH 3-150  
Without thermoblock  
Heating and cooling thermostat

Combitherm-2 CH 3-150 is specially designed to thermostabilise materials at temperatures from  $-3^{\circ}\text{C}$  to  $+150^{\circ}\text{C}$  according to methods of analysis. To obtain useful functionality and decrease foot-print of instruments Combitherm-2 ...

[read more](#)