

## B10-16, block

### DESCRIPTION

Ø16 mm x 10 sockets depth 56 mm

### SPECIFICATIONS

Depth	56 mm
Capacity	10
Diameter	16 mm



### CAT. NUMBER

BS-010418-BK

B10-16



CH 3-150  
Without thermoblock  
Heating and cooling thermostat

Combitherm-2 CH 3-150 is specially designed to thermostabilise materials at temperatures from  $-3^{\circ}\text{C}$  to  $+150^{\circ}\text{C}$  according to methods of analysis. To obtain useful functionality and decrease foot-print of instruments Combitherm-2 ...

[read more](#)