

B23-1.5, block

DESCRIPTION

23 sockets for 1.5 ml microtest tubes depth 35 mm

SPECIFICATIONS

Depth	35 mm
Amount of 1.5 ml tubes	23



CAT. NUMBER

BS-010418-DK

B23-1.5



CH 3-150
Without thermoblock
Heating and cooling thermostat

Combitherm-2 CH 3-150 is specially designed to thermostabilise materials at temperatures from -3°C to $+150^{\circ}\text{C}$ according to methods of analysis. To obtain useful functionality and decrease foot-print of instruments Combitherm-2 ...

[read more](#)