

## External temperature probe


### DESCRIPTION

An external probe provides direct control of the stirred liquids temperature.

The cable is covered with Teflon, mechanically strong, elastic and chemically stable against oils, acids, aggressive reagents and liquids

### SPECIFICATIONS

Operation temperature range	-50 °C to +250 °C
Cable length	1 m
Connection	type K
Probe type	Thermocouple



### CAT. NUMBER

BS-010309-BK

External temperature probe



Intelli-Stirrer MSH-300i  
BS-010309-AAA  
Magnetic Stirrer with hot plate

Intelli-Stirrer MSH-300i is magnetic stirrer of the new generation. Enclosures of stirrers are made of metal painted with powder enamel chemically resistant to acids and alkali. The stirrers are equipped ...

[read more](#)