

FC-2000, clamp

DESCRIPTION

Clamp 2000 ml - Ø165 mm for UP-168



CAT. NUMBER

BS-010126-NK

FC-2000, clamp 2000ml

ACCESSORIES



UP-168
BS-010135-JK
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



PSU-20i
Without platform
Multi-functional Orbital Shaker

Shaker PSU-20i provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)