


Glass test tubes 18mm

DESCRIPTION

Glass Test Tubes 18x100mm, high borosilicate, PP Cap with silicone pad.
Packing - 100 pcs/box

SPECIFICATIONS

Height	100 mm
Diameter	18 mm



CAT. NUMBER

BS-050102-NK Glass sample tubes with lids 18x100, 100 pcs.

ACCESSORIES



DEN-1
BS-050102-AAF
Densitometer (suspension turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units (0 – 180×10^7 cells/ml).

Densitometer provides the opportunity to measure solution ...

[read more](#)



DEN-1B
BS-050104-AAF
Densitometer (suspension turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units (0 – 180×10^7 cells/ml).

Densitometers provide the opportunity to measure solution turbidity in ...

[read more](#)