

LF-1, Drawer unit

DESCRIPTION

New modular design of laboratory furniture provides flexibility and ease of use.

SPECIFICATIONS

Drawer internal dimensions (W×D×H)	235 × 385 × 50 mm</br>235 × 385 × 150 mm
Overall dimensions (W×D×H)	300 × 450 × 705 mm
Weight	28 kg




CAT. NUMBER

BS-050101-BK

LF-1, Drawer unit

ACCESSORIES



UVC/T-AR

BS-040102-AAA

DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box UVC/T-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, plexiglas walls and ...

[read more](#)



UVC/T-M-AR

With built in socket and inlet for power cords
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box UVC/T-M-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)




LT-120

BS-040107-EK

Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



UVT-B-AR

With built in socket and inlet for power cords
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box UVT-B-AR is designed for clean operations with DNA samples. UV-Cabinet provide protection against contamination.

Model is a bench-top type, made of metal framework, working surface made of ...

[read more](#)



UVT-S-AR

BS-040107-AAA

DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box UVT-S-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



LT-150

BS-040107-JK

Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



LT-180

BS-040107-HK

Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



T-4

BS-040101-BK

Table

New modular design of laboratory furniture provides flexibility and ease of use.



T-4L

BS-040107-BK

Table

New modular design of laboratory furniture provides flexibility and ease of use.