

P-02/05, platform

DESCRIPTION

For 24 tubes 0.5 ml and 48 tubes 0.2 ml.

SPECIFICATIONS

Capacity

24 and 48

Microtubes volume

0.5 ml and 0.2 ml



CAT. NUMBER

BS-010216-DK

Platform P-02/05



MPS-1
With built-in universal platform
High-Speed Multi Plate Shaker

High-Speed Multi Plate Shaker **MPS-1** can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes ...

[read more](#)