

## P-02/96, platform

### DESCRIPTION

Platform for 96 tubes 0.2 ml or semi-/ unskirted PCR plate.

### SPECIFICATIONS

Capacity	96
Plate	Semi-/unskired PCR plate
Microtes tubes volume	0.2 ml



### CAT. NUMBER

BS-010216-CK

Platform P-02/96



**MPS-1**  
With built-in universal platform  
High-Speed Multi Plate Shaker

High-Speed Multi Plate Shaker **MPS-1** can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes ...

[read more](#)