

# P-6/1000, platform

## DESCRIPTION

Platform with clamps for flasks.

## SPECIFICATIONS

Working area	360 × 400 mm
Flask volume	1000 ml
Number of clamps	6



## CAT. NUMBER

BS-010135-DK

Platform P-6/1000

## ACCESSORIES



**ES-20/60**  
Without platform  
Orbital Shaker-Incubator

Orbital Shaker-Incubator ES-20/60 for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



**PSU-20i**  
Without platform  
Multi-functional Orbital Shaker

Shaker PSU-20i provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)