

PRS-8/22, platform

DESCRIPTION

8 tubes 20-30 mm diameter (50 ml tubes), 22 tubes 10-16 mm diameter (1.5 ml-15 ml tubes).

SPECIFICATIONS

Platform type	optional
Tube diameter	20-30 and 10-16 mm
Tube volume	up to 50 and 1.5-15 ml
Capacity	14



CAT. NUMBER

BS-010118-AK

Platform PRS-8/22



Multi RS-60
Including PRS-48 platform
Programmable rotator

Programmable Rotator **Multi RS-60** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as well ...

[read more](#)