

PV-32, platform

DESCRIPTION

V-32 is supplied with a 32-socket universal platform for Eppendorf type tubes up to 1.5 ml (1.5/0.5/0.2 ml — 16/8/8 sockets) PV-32 and PL-1 platform for mixing single tube up to 50 ml.

SPECIFICATIONS

Tube volume	1.5/0.5/0.2 ml
Capacity	16/8/8



CAT. NUMBER

BS-010207-CK

PV-32



V-32
Including two platforms PV-32, PL-1
Multi-Vortex

Multi-Vortex V-32 is intended for intensive stirring of bacterial and yeast cell, washing from the culture medium and extraction of metabolites and enzymes from cells and cell cultures.

It is ...

[read more](#)