

# PV-48, platform

## DESCRIPTION

Platform for 6-8 × 0.2 ml strips or 48 tubes of 0.2 ml.

## SPECIFICATIONS

Strips	0.2 ml
Capacity	6-8



## CAT. NUMBER

BS-010207-GK

PV-48



V-32  
Including two platforms PV-32, PL-1  
Multi-Vortex

Multi-Vortex V-32 is intended for intensive stirring of bacterial and yeast cell, washing from the culture medium and extraction of metabolites and enzymes from cells and cell cultures.

It is ...

[read more](#)