

R-0.5/0.2M, rotor

DESCRIPTION

Rotor for 12 × 0.5 ml and 12 × 0.2 ml microtest tubes

SPECIFICATIONS

Amount of 0.5 ml tubes	12
Amount of 0.2 ml tubes	12



CAT. NUMBER

BS-010201-BK

Rotor R-0.5/0.2M



FV-2400 Micro-Spin
Including rotors R-1.5M, R-0.5/0.2M
Mini-Centrifuge/Vortex

Mini-centrifuge/vortex Micro-Spin **FV-2400** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories.

Micro-Spin provides simultaneous mixing ...

[read more](#)