

# R-12/10, Rotor

## DESCRIPTION

Swing-out rotor for round bottom centrifuge tube.

## SPECIFICATIONS

Rotor type	Swing-out
Tube dimensions (Ø × length)	16 × 105 mm
Capacity	12
Tube's volume	10–15 ml
Max. speed	4200 rpm
Max. RCF: LMC-3000	1610 × g
Max. RCF: LMC-4200R	3160 × g
Tube manufacturers	Nunc, Greiner, Greiner Bio-one, TPP etc.
Autoclavable	+



## CAT. NUMBER

BS-010208-BK

Rotor R-12/10



**LMC-3000**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-12**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-4200R**  
Without rotor  
Laboratory Refrigerated  
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**BN-13/75**  
BS-010208-PK  
adapter set

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-13/75 for vacutainers 2-5 ml



**BN-13/100**  
BS-010208-QK  
adapter set

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-13/100 for vacutainers 4-8 ml



**BN-16/100**  
BS-010208-RK  
adapter set

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-16/100 for vacutainers 8-10 ml



**RR-U**  
BS-010208-UK  
rack for rotors

Rack for rotors



**BN-16/90**  
BS-010208-SK  
bucket set for rotor R-12/10 R-24/10

Spare bucket set for rotor R-12/10, R-24/10. Plastic buckets (12 pcs).