

# R-12/15, Rotor

## DESCRIPTION

Swing-out rotor for conical bottom centrifuge tube.

## SPECIFICATIONS

Rotor type	Swing-out
Tube dimensions (Ø × length)	17 × 120 mm
Capacity	12
Tube's volume	15 ml
Max. speed	4200 rpm
Max. RCF: LMC-3000	1610 × g
Max. RCF: LMC-4200R	3160 × g
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



## CAT. NUMBER

BS-010208-EK Rotor R-12/15

## ACCESSORIES



**LMC-3000**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-12**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-4200R**  
Without rotor  
Laboratory Refrigerated Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**RR-U**  
BS-010208-UK  
rack for rotors

Rack for rotors



**BN-17/120**  
BS-010208-TK  
bucket set for rotor R-12/15

Spare bucket set for rotor R-12/15. Plastic buckets (12 pcs).