

## R-1.5M, rotor

### DESCRIPTION

Rotor for 12 x 1.5/2 ml microtest tubes

### SPECIFICATIONS

Amount of 2 ml tubes	12
Amount of 1.5 ml tubes	12



### CAT. NUMBER

BS-010201-AK	Rotor R-1.5M
--------------	--------------



**FV-2400 Micro-Spin**  
Including rotors R-1.5M, R-0.5/0.2M  
Mini-Centrifuge/Vortex

Mini-centrifuge/vortex Micro-Spin **FV-2400** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories.

Micro-Spin provides simultaneous mixing ...

[read more](#)