

Ri-24/10, Rotor

DESCRIPTION

Swing-out rotor for round bottom centrifuge tube.

SPECIFICATIONS

Rotor type	Swing-out
Tube dimensions (Ø × length)	16 × 105 mm
Capacity	24
Tube's volume	10–15 ml
Max. speed	4000 rpm
Max. RCF	2860 × g
Tube manufacturers	Nunc, Greiner, Greiner Bio-one, TPP etc.
Autoclavable	+



CAT. NUMBER

BS-010221-LK

Rotor Ri-24/10



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BN-13/75-24
BS-010221-RK
adapter set

For rotors R-24/10, Ri-24/10 optional adapter set (24 pcs) - BN-13/75-24 for vacutainers 2-5 ml



BN-13/100-24
BS-010221-SK
adapter set

For rotors R-24/10, Ri-24/10 optional adapter set (24 pcs) - BN-13/100-24 for vacutainers 4-8 ml



BN-16/100-24
BS-010221-TK
adapter set

For rotors R-24/10, Ri-24/10 optional adapter set (24 pcs) - BN-16/100-24 for vacutainers 8-10 ml



RR-U
BS-010208-UK
rack for rotors

Rack for rotors



BN-16/90
BS-010208-SK
bucket set for rotor R-12/10 R-24/10

Spare bucket set for rotor R-12/10, R-24/10. Plastic buckets (12 pcs).