

# Ri-24GC, Rotor

## DESCRIPTION

Rotor for 24 Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen)

## SPECIFICATIONS

Rotor type	Swing-out
Card dimensions (l × w)	53 × 74 mm
Capacity	24 gel cards
Max. speed	1500 rpm
Max. RCF	280 × g
Gel card Manufacturers	Grifols, DiaMed, etc.
Autoclavable	+



## CAT. NUMBER

BS-010221-NK

Rotor Ri-24GC

## ACCESSORIES



**LMC-56**  
Without rotor  
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**RR-U**  
BS-010208-UK  
rack for rotors

Rack for rotors