

# Ri-6P, Rotor

## DESCRIPTION

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).

## SPECIFICATIONS

|                              |  |
|------------------------------|--|
| Tube holders material        | POM Kocetal  |
| Rotor type                   | Swing-out  |
| Tube dimensions (Ø × length) | 29 × 115 mm  |
| Capacity                     | 6  |
| Tube's volume                | 50 ml  |
| Max. speed                   | 4200 rpm   |
| Max. RCF                     | 3160 × g   |
| Max. temperature             | 150 °C   |
| Tube manufacturers           | Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc. |
| Autoclavable                 | +  |



## CAT. NUMBER

BS-010221-IK

Rotor Ri-6P

## ACCESSORIES



**LMC-56**  
 Without rotor  
 Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**BI-25-6**  
 BS-010221-VK  
 adapter set

For rotors R-6, R-6P, Ri-6, Ri-6P optional insert set (6 pcs) - BI-25-6 for 25 ml conical tubes



**RR-U**  
 BS-010208-UK  
 rack for rotors

Rack for rotors



**BN-11/30**  
 BS-010208-ZK  
 bucket set for rotor R-6P

Spare bucket set for rotor R-6P. Plastic buckets (POM Kocetal) (6 pcs).