

RMT-24, Rotor

DESCRIPTION

Fixed-angle rotor for 24 - 1.5 or 2 ml microtubes

SPECIFICATIONS

Capacity	24
Max. speed	6000 rpm
Max. RCF	3750 × g



CAT. NUMBER

BS-010221-BK

Rotor RMT-24

ACCESSORIES



LMC-4200R
Without rotor
Laboratory Refrigerated
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



RR-U
BS-010208-UK
rack for rotors

Rack for rotors