

RR-U, rack for rotors

DESCRIPTION

Rack for rotors



CAT. NUMBER

BS-010208-UK

Rack for rotors RR-U

ACCESSORIES



LMC-3000
Without rotor
Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-12
Without rotor
Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-4200R
Without rotor
Laboratory Refrigerated
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BR-4U
BS-010221-AK
Rotor

Bucket rotor without bucket adapters



BA-2/50
BS-010221-CK
Adapters

Bucket adapters for 2 tubes of 50 ml (set of 4 pcs.)



BA-4/15
BS-010221-EK
Adapters

Bucket adapters for 4 tubes of 15 ml (set of 4 pcs.)



BA-4/10
BS-010221-DK
Adapters

Bucket adapters for 4 tubes of 10 ml (set of 4 pcs.)



BA-8/5
BS-010221-FK
Adapters

Bucket adapters for 8 tubes of 5 ml (set of 4 pcs.)



BA-14/2U
BS-010221-JK
Adapters

Bucket adapters for 14 microtest tubes of 2 ml (set of 4 pcs.)



R-6
BS-010208-DK
rotor

Swing-out rotor for conical bottom centrifuge tube. Aluminium buckets.



R-6P
BS-010208-XK
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).



Ri-6
BS-010221-HK
rotor

Swing-out rotor for conical bottom centrifuge tube. Aluminium buckets.



Ri-6P
BS-010221-IK
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).



RMT-24
BS-010221-BK
Rotor

Fixed-angle rotor for 24 - 1.5 or 2 ml microtubes



R-12/10
BS-010208-BK
Rotor

Swing-out rotor for round bottom centrifuge tube.



Ri-12/10
BS-010221-MK
Rotor

Swing-out rotor for round bottom centrifuge tube.



R-24/10
BS-010212-JK
Rotor

Swing-out rotor for round bottom centrifuge tube.



Ri-24/10
BS-010221-LK
Rotor

Swing-out rotor for round bottom centrifuge tube.



R-12/15
BS-010208-EK
Rotor

Swing-out rotor for conical bottom centrifuge tube.



R-2
BS-010208-AK
Rotor

Swing-out rotor for standard 96-well microtitre plates, skirted PCR plates, deepwell plates up to 45mm.



Ri-12/15
BS-010221-KK
Rotor

Swing-out rotor for conical bottom centrifuge tube.



Ri-2
BS-010221-GK
Rotor

Swing-out rotor for standard 96-well microtitre plates, skirted PCR plates, deepwell plates up to 45mm.



R-24GC
BS-010208-VK
Rotor

Rotor for 24 Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen)



Ri-24GC
BS-010221-NK
Rotor

Rotor for 24 Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen)