

# SC-96AC, block

## DESCRIPTION

96-well unskirted or semi-skirted microplate (0.2 ml) for PCR or 12 × 8 - 0.2ml strips or 96 tubes of 0.2 ml.

## SPECIFICATIONS

96-well unskirted or semi-skirted microplate (0.2 ml)

1



## CAT. NUMBER

BS-010143-FK

SC-96AC

## ACCESSORIES



### TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker TS-100C provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



### TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker TS-100C Smart provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)