

SV-4/30, platform

DESCRIPTION

4 sockets for 50 ml (Ø 30 mm) tubes

SPECIFICATIONS

Diameter	30 mm
Tube volume	50 ml
Capacity	4



CAT. NUMBER

BS-010210-AK

SV-4/30



MSV-3500
With all 4 platforms included
Multi Speed Vortex

Multi Speed Vortex **MSV-3500** is designed for soft or intensive mixing of reagents in different size and type plastic tubes (0.2 to 50 ml).

It is designed for operation in ...

[read more](#)