

TR-30/13, rack

DESCRIPTION

Test tube rack for \varnothing 10 to \varnothing 13mm tubes, capacity 30 tubes. WxDxH: 155x90x112mm

SPECIFICATIONS

Capacity	30
Tube size	\varnothing 10-13 mm



CAT. NUMBER

BS-010406-1K

Test tube rack for \varnothing 10 to \varnothing 13mm tubes, capacity 30 tubes



WB-4MS
With base BP-1 and lid
Stirred water bath

Water bath-thermostat **WB-4MS** is designed for chemical, pharmaceutical, medical and biological laboratory research.

WB-4MS provides increased temperature stabilization (up to 0.1°C) due to built-in magnetic stirrer (regulated speed 250–1,000 rpm).

...

[read more](#)