

TR-5/30, rack

DESCRIPTION

Test tube rack for $\varnothing 30$ mm tubes, capacity 5 tubes. WxDxH: 155x90x112mm

SPECIFICATIONS

Capacity	5
Tube size	$\varnothing 30$ mm



CAT. NUMBER

BS-010406-KK	Test tube rack for $\varnothing 30$ mm tubes, capacity 5 tubes
---------------------	--



WB-4MS
With base BP-1 and lid
Stirred water bath

Water bath-thermostat **WB-4MS** is designed for chemical, pharmaceutical, medical and biological laboratory research.

WB-4MS provides increased temperature stabilization (up to 0.1°C) due to built-in magnetic stirrer (regulated speed 250–1,000 rpm).

...

[read more](#)