

UP-12, platform

DESCRIPTION

Universal platform with adjustable bars for different types of flasks, bottles, and beakers with silicone mat.

SPECIFICATIONS

Working area	270 × 185 × 40 mm
Flask volume	universal



CAT. NUMBER

BS-010108-AK

Platform UP-12

ACCESSORIES



HB-200
BS-010108-FK

Additional holding bar for UP-12



ES-20
Without platform
Orbital Shaker-Incubator

The **ES-20** is a compact benchtop Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program

Built-in ...

[read more](#)



PSU-10i
Without platform
Orbital Shaker

Shaker **PSU-10i** provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5) ...

[read more](#)



CPS-20
Without platform
CO₂ Shaker

CO₂ Shaker **CPS-20** provides regulated orbital motion of the platform and is designed for use specifically in CO₂ incubators. CPS-20 is specifically designed for use in harsh environments such as ...

[read more](#)